

# Practice Information:

## Supporting therapy decision in Hypertension




### Clinical Background

Often the question is raised if one should start intervention on a patient presenting a practice blood pressure of 130-159/85-99 mmHg or a LDL-value of 130 – 190 mg/dl (3.4 – 4.9 mmol/l).

The decision to start intervention can be based on evidence of subclinical end-organ-damage such as arterial stiffness measured by aortic pulse wave velocity (aPWV).

Evidence based data from Ben-Shlomo et al. demonstrate that aPWV is predictive for cardio-vascular events, independent of age. The predictive value of aPWV is significantly superior to the sole capture of common risk factors, such as hypertension, hypercholesteremia or the Framingham risk score.




**Fig. 1**  
Differentiation between high and low risk patients (in percentage)

### Therapy Management

**PWV measurement:** The use of a standard upper-arm cuff allows the simultaneous measurement of blood pressure and aPWV in one process stratification: An increased aPWV value classifies the patient into a “high risk”, instead of “low risk” group.<sup>2</sup>

**Clinical Monitoring** Therapy progress and success can be documented by the measurement of aPWV.



**Fig. 2**  
aPWV re-classifies patients according to their individual risk

### Reimbursement

aPWV measurement may be offered to patients as a “vascular age measurement”.

The measurement can be utilized within the framework of prevention programs (for example during check-up programs).

Reimbursement for aPWV measurement is expanding to many countries. In others, the examination is commonly paid for by the patients themselves.

# Practice Information:


## Supporting therapy decision in Hypertension

### Methodology

**Device Technology:** Mobil-O-Graph® (Mobil-O-Graph® NG Classic including PWA Upgrade or Mobil-O-Graph® 24h PWA).

**Measurement:** Blood pressure and aortic PWV measurement in one process.

**Clinical Validation:** The measurement of aPWV using Mobil-O-Graph® correlates significantly with invasive catheter measurements. The very high correlation of  $R = 0.81$  is also confirmed by a reproducibility of 0.05 m/s (ranging from -0.47 to 0.57 m/s).<sup>3</sup>



**Fig. 3**  
High agreement between Mobil-O-Graph® and invasive catheter measurement of PWV

### Technology Features

**Measurements:** To be applied as Practice/Clinic or as 24h ambulatory measurement.

**HMS-CS Upgrade:** The HMS CS software application can be easily upgraded to perform practice/clinic pulse wave velocity by using the licensed PWA Dongle (USB).



**Fig. 4 Upgrade-Dongle**  
HMS CS software upgrade PWA Dongle in combination with Mobil-O-Graph® NG Classic 24h ABPM enables the measurement of blood pressure and pulse wave velocity in one process

### Further Information

For more detailed information on the measurement of aPWV measurement in your practice and/or clinic we invite you to visit our I.E.M. GmbH web-site ([www.iem.de](http://www.iem.de)) or request information under [info@iem.de](mailto:info@iem.de) or via + 49 2402 9500-0.

#### Literature and sources:

1. Aortic Pulse Wave Velocity Improves Cardiovascular Event Prediction; An Individual Participant Meta-Analysis of Prospective Observational Data From 17,635 Subjects; Yoav Ben-Shlomo, et al ; Journal of the American College of Cardiology Vol. 63, No. 7, 2014
2. 2013 ESH/ESC Guidelines for the management of arterial hypertension; Journal of Hypertension 2013, 31:1281–1357 (see section 4.2.3)
3. Oscillometric estimation of aortic pulse wave velocity: comparison with intra aortic catheter measurements; Weber T, Hametner B, et al; Blood Press Monit. 2013 Jun;18(3):173-6.