Mobil-O-Graph® PWA

Gebrauchsanweisung	DE
Instructions for Use	EN
Návod k použití	CS
Brugsanvisning	DA
Εγχειρίδιο Λειτουργίας	EL
Manual de instrucciones	ES
Käyttöohje	FI
Mode d’emploi	FR
Upute za uporabu	HR
Istruzioni per l'uso	IT
작동 매뉴얼	KO
Naudojimo instrukcijos	LT
Gebruiksaanwijzing	NL
Bruksanvisning Instrukcja	NO
obsługi Instrucțiuni	PL
Instruções de utilização	PT
Instrucțiuni de utilizare	RO
utilizare Bruksanvisning	SV
Kullanım kılavuzu	TR

IEM GmbH
Gewerbepark Brand 42
52078 Aachen
Germany

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Instructions for Use
7.1_2026-09-14
EN

Mobil-O-Graph® PWA

Long-term blood pressure monitor and pulse wave analysis monitor (PWA)

These instructions for use describe the following device variants (GTIN): 04041346102269, 04041346102337, 04041346000527

Note

For USA: Federal law restricts the sale of this device to be by or on the order of a doctor.

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The contents of this technical manual must not be duplicated or published without the written authorization of IEM GmbH.

© IEM GmbH 2026. All rights reserved.

Table of Contents

1 Introduction ... 5

1.1 About these instructions for use ... 5

1.2 Clinical trials ... 5

1.3 CE marking ... 5

1.3.1 Applicable standards and directives ... 5

2 Instruction notes ... 7

2.1 Intended purpose ... 7

2.2 Essential performance features ... 7

2.3 Indication and contraindication ... 7

2.3.1 Indication ... 7

2.3.2 Contraindication ... 8

2.4 Side effects of long-term blood pressure measurement ... 8

2.5 Safety ... 8

2.5.1 Defining the signal words used ... 8

2.6 Important safety instructions for the doctor ... 9

2.7 Important safety instructions for the device ... 14

3 Product description ... 16

3.1 Description ... 16

3.2 Unpacking ... 16

3.2.1 Approved accessories ... 16

3.3 Description of device ... 17

3.3.1 The buttons ... 17

3.3.2 The display ... 18

3.3.3 The acoustic signals ... 18

3.3.4 Cuff connection ... 19

3.3.5 The data socket ... 19

3.4 Preparation for measurement ... 19

3.4.1 Switching on ... 19

3.4.2 Clearing the memory ... 19

3.4.3 Set time/date ... 20

3.4.4 Transfer of patient data ... 20

3.4.5 Putting on the blood pressure monitor and starting the measurement ... 20

3.5 Positioning and behavior of the patient ... 23

3.6 Technical data and ambient conditions ... 23

3.7 Symbols ... 24

4 Care and maintenance ... 26

4.1 Cleaning ... 26

4.2 Disinfection ... 28

4.3 Maintenance schedule ... 29

5 Troubleshooting ... 30

5.1 Main sources of error....30
5.2 Error messages....31
6 Rechargeable batteries/batteries....32
6.1 Operation with rechargeable batteries....32
6.2 Operation with ALKALINE batteries....32
7 Disposal....33
8 Warranty and repair conditions....34
Appendix....35
A Appendix 1: EMC Guidelines and Manufacturer's Declaration....36

1 Introduction

1.1 About these instructions for use

These instructions for use explain the blood pressure monitor and the accessories in the order in which you will put the device into operation and use it later.

The individual functions will only be explained when you actually need them. You will therefore be familiarized with the blood pressure monitor step by step.

The instructions for use must be kept for later use and must always be available!

1.2 Clinical trials

The Mobil-O-Graph® PWA meets the requirements of the ESH (European Society of Hypertension), BSH (British Hypertension Society) and ISO 81060-2:2018+AMD1:2020+AMD2:2024.

1.3 CE marking

CE 2797

The Mobil-O-Graph® PWA fulfills the requirements of the following standards and directives (1.3.1 Applicable standards and directives, p. 5) and bears the CE marking.

IEM GmbH hereby declares that the Mobil-O-Graph® PWA complies with the labeling of a medical device according to Regulation (EU) 2017/745.

Tip

The full text of the EU Declaration of Conformity is available at the following internet address: www.iem.de/doc/

As part of product monitoring, IEM GmbH is interested in information and experiences gathered from using the blood pressure monitor and the instruction manual. We are therefore grateful for any relevant feedback.
If you are unsure about anything set out in this instruction manual, IEM GmbH will be happy to answer any queries.
Serious incidents occurring in connection with the product must be reported to IEM GmbH and, where appropriate, the responsible authority of the respective member state immediately.

1.3.1 Applicable standards and directives

- Regulation (EU) 2017/745 (MDR)
- Regulation (EU) 2021/2226
- Regulation (EU) 2025/1234
- Directive 2011/65/EU (RoHS)
- IEC 60601-1: 2005+AMD1:2012+AMD2:2020
- IEC 60601-1-2: 2014+AMD1:2020
- EC 60601-1-6: 2010+AMD1:2013+AMD2:2020
- IEC 60601-1-11:2015+AMD1:2020

- IEC 80601-2-30:2018
- ISO 81060-2:2018+AMD1:2020+AMD2:2024
- FDA 510(K)

2 Instruction notes

2.1 Intended purpose

The intended use of the Mobil-O-Graph® PWA is 24-hour blood pressure measurement. Users have the option of activating a PWA license so that the device also records pulse wave data.

Note

The Mobil-O-Graph® PWA should only be used under the medical supervision of qualified medical staff.

Note

There are currently no clinical studies concerning the use of pulse wave analysis in children. Accordingly, there is not yet a confidence interval for people under the age of 20.

2.2 Essential performance features

The essential performance features are defined as blood pressure measurement with:

- Error tolerances of the pressure gauge and measurement results within the required limits according to IEC 80601-2-30
- Maximum change value in blood pressure determination according to IEC 80601-2-30
- Power delivery (pressure supply to the cuff) within the set limits according to IEC 80601-2-30
- An error is issued if successful blood pressure measurement is not possible.

The device does not trigger any alarms in the sense of IEC 60601-1-8 and is not intended for use together with HF surgical devices or for clinically monitoring patients, e.g. in an intensive care unit.

In the event that the status or configuration of the device is unclear, the device will go into a safe operating mode, which causes the device to release the air from the cuff. The cuff is not automatically pressurised; to do so, the device must be initiated manually.

2.3 Indication and contraindication

2.3.1 Indication

The Mobil-O-Graph® PWA is intended to clarify the blood pressure situation and for diagnosis support. The field of application is home healthcare environment and professional healthcare facility environment such as doctor's practices and medical centres. The doctor may, at his own discretion, carry out an appropriate examination on his patients with this medical system if, among other things, they:

- are hypotensive or
- are hypertensive,
- require antihypertensive therapy,
- suffer from myocardial hypertrophy or
- suffer from nephrological dysfunction.

2.3.2 Contraindication

- The Mobil-O-Graph® PWA is not intended for use with newborn babies and children under the age of 3!
- The Mobil-O-Graph® PWA is not intended for independent use by people with cognitive impairments and unsupervised children!
- The cuffs are not suitable for patients who are hypersensitive (allergic reactions) to epoxy resin!
- The Mobil-O-Graph® PWA is not intended for use on limbs with an intravascular access or intravascular treatment, or with an arteriovenous (AV) shunt!
- The Mobil-O-Graph® PWA is not intended for use on limbs on the same side where a mastectomy has been performed!
- The Mobil-O-Graph® PWA with PWA functionality is not intended for use in emergency situations!

2.4 Side effects of long-term blood pressure measurement

Long-term blood pressure measurement is a commonly practised and valued measurement technique,and has found its way into daily diagnostics and treatment monitoring.

Blood pressure measurement can lead to petechial bleeding or haematomas on the measuring arm despite the correct positioning of the cuff in case of coagulation disorders, taking anticoagulants or in the case of sensitive body tissue. Always check whether the patient has coagulation disorders or is undergoing anticoagulant treatment. The patient-dependent risk resulting from anticoagulant treatment or in patients with coagulation disorders arises irrespective of the type of measuring device.

2.5 Safety

Read the safety instructions carefully before using the products! It is important that you understand the information in these instructions for use. If you have any questions, please do not hesitate to contact technical support.

2.5.1 Defining the signal words used

The following symbols and signal words are used in these instructions for use to draw attention to dangers and important information:

⚠ DANGER

Highest category of safety instructions!

The severity of damage ranges from moderate to serious or catastrophic. Failure to observe the safety instructions will result in the following:

• serious injury or death of the patient, user or third parties, or

• serious environmental damage.

This warning symbol in conjunction with the signal word **DANGER** indicates an imminent hazard with a frequent/occasional probability of occurrence.

⚠ WARNING

Second-highest category safety information!

The severity of the damage/injury ranges from moderate to serious, through to catastrophic. Failure to observe information may result in

• severe injuries or death,

• permanent disability/injury; permanent disability/injury in spite of medical intervention.

This warning symbol combined with the signal word **WARNING** indicates a possible imminent risk with rare/unlikely probability of occurrence.

⚠ CAUTION

Medium category of safety instructions!

The severity of damage ranges from moderate to minor. The following may result:

• minor or moderately serious injuries;

• short-term invalidity, possibly requiring uncomplicated and minor medical intervention, with no lasting invalidity;

• harm requiring complex medical intervention;

• minor misdiagnoses with negligible implications for medication or misdiagnoses with tangible implications for medication (dizziness, hyperactivity).

This warning symbol in conjunction with the signal word **CAUTION** indicates a potential hazard with a frequent/occasional/rare/unlikely probability of occurrence.

NOTICE

Short description of the danger!

Non-adherence may result in material damage to the product, its accessories or other property.

This warning symbol in conjunction with the signal word **NOTICE** indicates possible material damage.

Note

The signal word **Note** indicates medical information relating to the product.

Tip

The signal word **Tip** indicates:

- further information on the product or its accessories.
- references to external documents in which further information can be optionally found.

2.6 Important safety instructions for the doctor

⚠ DANGER

Risk of strangulation due to shoulder strap and cuff tube!

If handled incorrectly, the shoulder strap and cuff tube may wrap around the neck, impede breathing and result in severe injuries or even death.

• The shoulder strap and cuff tube must not be placed around the neck.

• The cuff must only be worn on the upper arm.

⚠ DANGER

• People (including children) who, due to their physical, sensory, or mental abilities, or their inexperience or lack of knowledge, are not able to use the device safely must not use the device without being supervised by a responsible person.

• For children, the device must be used with particular care and with constant supervision as prescribed by the doctor.

• People who are not responsible for their actions must not use the device. The device must be stored so that they cannot access it.

• In the event of pain, swelling, redness, or numbness during the measurement: switch off the device immediately, remove the cuff and inform your doctor.

• The measuring process can be stopped at any time by pressing any button. The cuff will be deflated and the device can be removed.

⚠ WARNING

Lymph congestion due to inflation of the cuff!

Inflating the cuff after lymph node removal increases the risk of lymph congestion. If there is existing lymphedema, the existing swelling can increase.

Only use the blood pressure monitor when taking the possible consequences into account.

⚠ WARNING

Allergic reactions due to the cuff material!

In very rare cases, the epoxy resin on the cuff material may cause allergic reactions in hypersensitive people, for example, skin redness and weeping pustules on the hands, legs, face, and neck.

• If you experience pain or allergic reactions, immediately remove the cuff.

• When positioning the blood pressure cuff, ensure that no disinfectant residue remains on it.

• Pay attention to hygiene and observe the information on cleaning and disinfection.

⚠ WARNING

Risk of injury due to incorrectly connected air tube!

An incorrectly connected air tube may impede the operation of the device and result in swelling or numbness on the arm around which the cuff is placed.

Do not bend, knot or stretch the air tube.

⚠ WARNING

Increased electromagnetic emissions/reduced electromagnetic immunity of the device due to the use of non-approved accessories!

The use of accessories, transducers and cables, other than those specified or supplied by the manufacturer of this device, may result in increased electromagnetic emissions or decreased electromagnetic immunity of this device and result in improper operation.

Only use approved accessories to ensure proper operation.

⚠ WARNING

Misdiagnosis due to incorrect measurements in patients with white coat hypertension!

In patients with white coat hypertension, stress and anxiety can lead to misdiagnosis due to false readings.

⚠ WARNING

Switching off the screen via the analysis unit can help patients with white coat hypertension to minimize stress and anxiety during 24-hour blood pressure measurement.

⚠ WARNING

Self-diagnosis by the patient

The patient derives therapeutic measures from the displayed measurement results.

Deactivate the display of the measurement results.

⚠ WARNING

Incorrect or inaccurate measurement results due to an incorrectly positioned cuff

Incorrect measured values can lead to a misdiagnosis.

Select the correct cuff size

Ensure the correct alignment of the cuff during a measurement

Make sure that the cuff is adequately tight.

The cuff should be fitted by trained personnel.

Inform the patient about the correct fit of the cuff and advise them to check the fit of the cuff regularly, but don't adjust it during a measurement.

⚠ CAUTION

Circulation problems due to constant cuff pressure or due to measurements performed too often

Incorrect positioning of the cuff or performing measurements too often may result in swelling or numbness on the arm around which the cuff is placed.

• Ensure that the cuff tube is positioned correctly and ensure that the cuff tube is not knotted, compressed, kinked, or separated.

• After a measurement, wait **at least 3 minutes** before starting a new measurement.

• If you experience pain, swelling, redness, or numbness on the arm around which the cuff is placed, inform your doctor.
(It should be expected that light to moderate discomfort may be experienced during a blood pressure measurement.)

• The measurement process can be stopped at any time by pressing any button. This will deflate the cuff and the cuff can be removed.

⚠ CAUTION

Danger of a temporary loss of function of a present electrical medical device as a result of putting on or pumping up a cuff if you are wearing a further electrical medical device for monitoring on the same limb.

Other medical electrical devices located on the same limb as the cuff may fail temporarily. This means that body functions can no longer be monitored.

Only apply the cuff if you are not wearing another medical electrical device on this arm.

⚠ CAUTION

Risk of injury when applying and inflating a cuff over a wound!

The pressure created when inflating the cuff can aggravate minor injuries or open wounds.

⚠ CAUTION

Do not place the Mobil-O-Graph® PWA over a wound or bandage.

⚠ CAUTION

Danger of blood flow disruptions as a result of putting on and pumping up a cuff on limbs with an intravascular drip or intra-vascular treatment or with an arteriovenous (AV) shunt.

If there is an intra-vascular access, an infusion (e.g. intravenous drip or blood transfusion) or an arteriovenous (AV) shunt on the same limb as the cuff, circulatory disorders may occur, which can lead to minor or moderate injuries.

Do not apply the cuff to the arm with an intra-vascular access, intravascular treatment or with an arteriovenous (AV) shunt. The other arm can be used for blood pressure measurement.

⚠ CAUTION

Risk of injury due to internal bleeding in the tissue or hematoma!

People may suffer impaired blood circulation, which may result in minor to moderate injuries to the upper arm or crook of the arm.

- Position the blood pressure cuff correctly.
- Only use the original blood pressure cuff.
- If you are taking blood clotting medication or you suffer from clotting problems, inform your doctor before using the blood pressure monitor.
- People with sensitive body tissue may experience bleeding into the tissue or hematoma despite using the cuff correctly and positioning it properly.

⚠ CAUTION

Risk of injury due to positioning and inflating a cuff on an arm on the side of the body where there has been a mastectomy!

Inflating the cuff may cause pain, trauma and/or other injuries to the arm on the side of the body where a mastectomy (breast removal) has been performed.

- Do not place the cuff on the arm on the same side of the body where a mastectomy has been performed.
- If the mastectomy has only been performed on one side of the body, use the arm on the other side of the body for blood pressure measurement.

NOTICE

Incorrect or inaccurate measurement results in pregnant women or patients with pre-eclampsia!

Due to stiff arteries during pregnancy and pre-eclampsia, correct measurement results cannot be guaranteed.

Only use the blood pressure monitor when taking the possible consequences into account.

NOTICE

Incorrect or inaccurate measurement results in patients with cardiac arrhythmia!

Correct measurement results cannot be guaranteed due to hardened arteries.

Only use the blood pressure monitor on patients with cardiac arrhythmia, e.g. caused by premature atrial or ventricular beats, atrial fibrillation, arteriosclerosis or sclerosis, taking into account the possible consequences.

NOTICE

Incorrect measurement results or damage to property due to using unapproved accessories!

• Incorrect results may be issued.

• The device may be rendered unusable.

Only use the original blood pressure cuff.

NOTICE

Incorrect measurement results due to movement artifacts!

The patient should relax as much as possible during the measurement and must not speak, except to communicate their discomfort.

Before the first measured value is recorded, 5 minutes must pass in silence.

For a 24-hour measurement, the patient should adopt one of the three positions (sitting, standing, lying down) during a blood pressure measurement if possible.

The measurements can be influenced by the measurement location, the position of the patient, exertion or the physiological condition of the person having the measurements taken.

Note

Patients with weak circulation (e.g. due to existing heart problems/weakness) may suffer from nausea during the measurement. The decision as to whether a 24-hour measurement should be performed depends on the individual case and is to be decided by the doctor in charge of treatment.

Note

Blood pressure measurement for patients with a left ventricular assist device (LVAD) is possible, but IEM does not guarantee the accuracy of the measurement.

Note

Damage to the device

- No liquid may enter the device. If you suspect that liquid has entered the device during cleaning or use, the device must no longer be used.
- If the device has been exposed to moisture, switch off the device and remove the batteries. It is essential that you inform technical support or return the device to your specialist supplier or the manufacturer.
- The device must not be operated in the vicinity of an MRI scanner or in the immediate vicinity of other medical electrical equipment.
- The device is not suitable for simultaneous use with high-frequency surgical devices.
- The device must not be electrically connected to a PC or other device if it is still attached to a patient.
- The device must not be used in airplanes.
- Do not drop the device and do not weigh it down with objects.
 If you see external damage or suspect damage, return the device to your specialist supplier or the manufacturer.
- Do not use the device directly next to other devices or stacked with other devices, as this could result in incorrect operation. If it still needs to be used in the above manner, this device and the other devices should be observed to ensure that they are working properly.

- The cuffs and hose are made of electrically non-conductive material. This protects the device against the effects of a defibrillator discharge. During a defibrillator discharge, the device itself must not touch the patient, as the device can be damaged by such a discharge and could result in incorrect values being displayed.
- The device meets all the requirements of the EMC standards, but it should not be exposed to strong electromagnetic fields, as this can lead to malfunctions outside the limit values. As a result, please note that the device should be at least 30 cm (12 inches) away from portable RF communication devices.

Note

The pulse wave analysis provides additional indicators for possible risks, but is not a sufficient indicator for individual diseases or therapy recommendations.

Note

Electrostatic discharges, such as discharges from synthetic textile clothing, can trigger a reboot of the device. The device starts in the last used operating status.

Note

The blood pressure monitor may only be used by the person who was given the blood pressure monitor by a doctor/healthcare provider in order to identify cardiological abnormalities.
The device must not be passed on to third parties because the transmitted measurements may be misinterpreted and result in incorrect treatment.

Tip

Would you like to support your patients with 24-hour blood pressure measurement and provide them with important information on using the device? You can find the patient information in DIN A4 format at www.iem.de/mobil-o-graph.

You can print out this patient information and offer it as a guide. Inform your patients in advance about the possible dangers when using the device to ensure safe operation.

2.7 Important safety instructions for the device

Housing

NOTICE

Loss of warranty!

Opening the housing voids any warranty.

Do not open the device.

Batteries/rechargeable batteries

⚠ WARNING

Damage to the device and property damage due to incorrect battery handling!

Incorrect, improper handling of the batteries/rechargeable batteries may result in the device catching fire. People in the immediate vicinity of the device may suffer burns, which may result in serious injuries or even death.

- Batteries/rechargeable batteries must not be thrown into the fire.
- Batteries/rechargeable batteries must never be subjected to high temperatures.
- Do not attempt to charge batteries.
- Batteries/rechargeable batteries must not be opened.

⚠ WARNING

• Batteries/rechargeable batteries must not be short-circuited.

⚠ CAUTION

Skin irritation due to leaking batteries/rechargeable batteries!

Fluid that leaks from the batteries due to incorrect handling may result in skin irritation.

• If you come into contact with the fluid, rinse it off with plenty of water.

• If the fluid gets into your eyes, do not rub your eyes; instead rinse them immediately with water for 10 minutes and consult a doctor without delay.

NOTICE

Leaking (rechargeable) batteries can cause damage to the device and property!

If the (rechargeable) batteries are stored for too long without the blood pressure monitor being used, this can cause the (rechargeable) batteries to leak.

Remove the (rechargeable) batteries from the battery compartment as soon as they have discharged and if the blood pressure monitor is not likely to be used for an extended period of time.

Tip

When charging rechargeable batteries, follow the instructions for use provided by the manufacturer of the used charger.

3 Product description

3.1 Description

The ambulatory blood pressure measurement system consists of two main components:

- the Mobil-O-Graph® PWA with various cuffs and other accessories.
- the analysis unit for analysis of the measurement results by the doctor.

The blood pressure monitor can be read out using the analysis unit. With the software, stored measurement results can be read into a PC, displayed on the screen in the form of graphics, lists and statistics and printed out.

The Mobil-O-Graph® PWA can be prepared for the next patient immediately afterwards. With a little practice, this procedure can be completed in just a few minutes. This allows the doctor to use the device around the clock every working day.

The Mobil-O-Graph® PWA with the analysis unit is designed to allow documentation and visualization of a blood pressure profile during the day and night. Additional parameters such as nocturnal values and blood pressure fluctuations are detected. This allows the doctor to prescribe an individual and optimal medical therapy and to monitor its success.

Tip

Please refer to the evaluation unit instructions for use for software operating instructions.

3.2 Unpacking

All parts included in the scope of delivery have been properly packed and checked for completeness and functionality. Should the product be incomplete or damaged, please inform the supplier immediately.

3.2.1 Approved accessories

Scope of delivery:

- Mobil-O-Graph® shoulder bag
- Connection cable (USB-C)
- IEM ABPM cuff (size S)
- IEM ABPM cuff (size M)
- IEM ABPM cuff (size L)
- 2 Akaline batteries (AA)
- Instructions for downloading the evaluation unit

Accessories not included in the scope of delivery:

- IEM ABPM cuff (size XS)
- IEM ABPM cuff (size XL)

3.3 Description of device

3.3.1 The buttons

All buttons are located on the front of the blood pressure monitor.

Note

You can cancel the measurement process prematurely with any button; the pressure in the cuff is then quickly released.

ON/OFF 	The blood pressure monitor is switched on and off using the ON/OFF button. No recorded data are lost when the device is switched off. After pressing the ON/OFF button to switch on, the device performs an initial self-test: • Display of the current battery voltage for **3 seconds** • Display test (display of all segments with countdown) • After successful completion, the time is displayed. The device is ready for operation. If the interface cable is connected before the blood pressure monitor is switched on, the following appears on the display **co** (= communication).
DAY/ NIGHT 	The DAY/NIGHT button can be used to separate the waking and sleeping phases during recording, which is important for statistics and graphical displays. Both the standard values and the frequency of the measurements are adapted to the day and night rhythm. ***Tip*** The exact specification regarding the effects on the printout can be found in the operating instructions for the analysis unit. If the day/night interval is permanently defined in the analysis unit, the patient can start an interval manually up to two hours before it begins by pressing the DAY/NIGHT button. The display shows **-d--d-.** Short: The patient is instructed to press the DAY/NIGHT button when they go to bed and again when they get up in the morning. This allows the measurement interval to be individually adapted to the patient and helps you to analyze the blood pressure profile. In addition to the interval adjustment, you will find corresponding notes on the printout. If this button is not pressed, the interval change takes place according to the set protocol.
EVENT 	The EVENT button is used to set markers during a measurement cycle. The display shows **-E--E-**. An additional measurement is triggered to record an event that influences the blood pressure. The patient should note the circumstances of the event, such as medication taken, in a personal measurement log so that they can discuss the events with you. Pressing the button marks the additional measurement via an internal code so that it is highlighted when read out with the analysis unit.
START 	The START button is used to trigger the 24-hour measurement. If the START button is pressed, the number of measurements registered so far appears on the display and the initial measurement is started. It differs from the automatic measurement according to the measurement protocol in that the cuff is inflated in stages. The required pressure of the cuff is determined, at which the systolic blood pressure value can be measured. This maximum required inflation pressure is saved and immediately "approached" by direct inflation during subsequent automatic measurements. **⚠ WARNING** **Incorrect measured values on first measurement**

	The values of the first measurement should be checked by the doctor for plausibility so that correct automatic measurements can be taken subsequently and the correct fit of the cuff is ensured. The instructions in the chapters 3.4 Preparation for measurement, p. 19 and 5 Troubleshooting, p. 30 should be followed in the event of an incorrect measurement.

3.3.2 The display

The LCD display is located on the housing top of the blood pressure monitor. It displays useful information for the doctor and the patient regarding measurement data, blood pressure monitor settings and measurement errors.

3.3.3 The acoustic signals

If an event is started, feedback is provided by acoustic signals. The acoustic signals used consist of single or multiple pulse sequences.

The following pulse sequences are output:

Number of sound pulses	To be heard at
1	Press the button: • ON/OFF (switching on and off) • START Sound pulse at the number of measurements on the display, start of measurement Sound pulse at end of measurement, display of measurement result -- **> except in the night interval** • START for at least **5 seconds** Sound pulse for the number of measurements on the display Sound pulse for **clr** (clear event) on the display EVENT for at least **2 seconds** Sound pulse after deleting all measurements • ON/OFF or START or DAY/NIGHT or EVENT Sound pulse for **-StoP-** (manual measurement abort) on the display **-->** Abort the measurement and release the pressure in the cuff • DAY/NIGHT Sound pulse for **-d--d-** (setting markers, sleep and wake phase) on the display • EVENT Sound pulse for **-E--E-** (setting markers during a measurement cycle, additional measurement) on the display • Plugging in the interface cable • Removing the interface cable • Measurement error
5	• System error (e.g. measurement abort)
Permanent	• Serious system errors (e.g. the cuff pressure is greater than 15 mmHg for more than 10 seconds outside the measurement)

3.3.4 Cuff connection

The cuff connection is located on the housing top of the Mobil-O-Graph® PWA. This metal plug is needed to connect the blood pressure monitor to the cuff via the cuff tube and the metal sleeve of the cuff.

Note

The metal plug (air coupling) must always engage with an audible "click". Otherwise there will be a leaky connection between the Mobil-O-Graph® PWA and tube which will lead to measurement errors.

3.3.5 The data socket

The data socket is located on the left-hand side of the housing. The supplied connection cable is connected to this socket. This is a USB-C connection.

Note
Older versions of the Mobil-O-Graph® PWA may be equipped with a 4-pin data socket. For further information about these device versions, please contact Customer Service.

3.4 Preparation for measurement

3.4.1 Switching on

Always check the condition of your blood pressure monitor first before giving it to a patient. This is done by observing the first notifications on the display of the blood pressure monitor shortly after switching it on.

The following should be displayed in this sequence:

Test	Display	Comment
Battery status (volts)	2.85	at least 2.6 V for NiMH rechargeable batteries at least 3.10 V for alkaline batteries
Display Segment Test	999:999 to 000:000	First the digits (999:999 to 000:000), then all other symbols of the LCD appear in succession. Check whether all segments are correctly and fully displayed (the complete program code is checked for correctness in the background).
Current 24-hour period	21:45	Between 00:00 and 23:59

If an error occurs during the internal test, the blood pressure monitor will show E004 on the display and an audible signal will sound. For safety reasons, the blood pressure monitor will be out of operation. Send the blood pressure monitor immediately to your specialist supplier or directly to IEM GmbH for repair.

3.4.2 Clearing the memory

The memory must be empty before every measurement, i.e. there must be no blood pressure data from the previous patient in the memory. If there are still values in there, these can be cleared using the delete function of the evaluation unit.

You can clear the memory manually by holding down the start button for more than 5 seconds. As it is held down, all segments of the LCD will be displayed first of all, then 1 tone burst will be emitted, the number of saved measurements will be displayed briefly and then "clr" will be displayed. If the event button is now pressed for more than 2 seconds within 5 seconds, all measurements will be deleted.

3.4.3 Set time/date

The Mobil-O-Graph® PWA has an internal backup cell that keeps the time running for **a minimum of 24 hours** even after the rechargeable batteries or batteries have been removed from the battery compartment. Check the time and date before each series of measurements.

Tip

The time and date can be set with the respective analysis unit.

⚠ CAUTION

Incorrect assignment of measurement results due to incorrect time and date settings!

Incorrect time and date settings can lead to incorrect assignment of the measurement results and possibly to misdiagnosis with incorrect medication

Set the correct time and date before starting a series of measurements so that the measurement results can be correctly assigned during assessment.

Ensure that there is sufficient voltage to guarantee operation of the blood pressure monitor.

Further information on the device that must be observed before passing it on to patients can be found in chapter 3.4.1 Switching on, p. 19.

Further information on rechargeable batteries/batteries can be found in chapter 6 Rechargeable batteries/batteries, p. 32.

3.4.4 Transfer of patient data

The blood pressure monitor is parameterized by transferring patient data via the evaluation unit, ensuring that the data can be correctly assigned during retrieval.

3.4.5 Putting on the blood pressure monitor and starting the measurement

Always place the device bag over your clothing. By varying the length of the belt, you can use it either as a hip belt or as a shoulder belt. Alternatively, a normal belt that matches the clothing can also be used. Now put on the cuff. The correct fit of the cuff is important for correct measurement.

We recommend putting the cuff on a bare upper arm. The cuff can also be worn over a thin shirt or blouse.




Fig. 1: Connecting the air hose to the device

Please note the following points:

1. The monitor must <u>not</u> be connected to other external devices when the monitor is connected!
2. The cuff must be positioned so that the pressure hose cannot be kinked at any point! The course of the pressure hose should ensure free movement of the upper arm and be guided over the neck to the other side of the body. The hose connection on the cuff must therefore be directed upwards.
3. Make absolutely sure that the artery symbol is on the brachial artery!
4. The lower edge of the cuff should be about 2 cm above the crook of the elbow!

5. The cuff should be placed on the upper arm. A simple test can be performed to check that it is fitting correctly: It should be possible to insert one or two fingers under the cuff.
6. The correct size of the cuff is important for correct blood pressure measurement. To ensure that reproducible values can be measured, standardized measurement conditions should be available, i.e. the cuff size should be adapted to the patient. The measuring tape included in the scope of delivery is used to measure the arm circumference in the middle of the upper arm and select the cuff:

Upper arm circumference	Cuff
14 – 20 cm (5.5 – 7.9 inch)	XS
20 – 24 cm (7.9 – 9.5 inch)	S
24 – 32 cm (9.5 – 2.6 inch)	M
32 – 38 cm (12.6 – 15 inch)	L
38 – 55 cm (15 – 21.7 inch)	XL

7. Connect the pressure hose of the cuff to the blood pressure monitor. Push the hose firmly onto the connection, the pressure hose must engage audibly; when releasing, simply pull back the outer metal ring of the plug.
8. The blood pressure monitor is now correctly applied and ready for measurement.

Once all the previous steps have been completed, the blood pressure monitor can be put into operation. First, a manual measurement is taken by pressing the "START" button. You can use this measurement to determine whether the blood pressure monitor is working properly.

If errors occur, check the correct procedure for setting up and installing the monitor and accessories again. If this does not help, repeat the commissioning procedure.

Only after a successful manual measurement can the patient be discharged until the device is removed.

Note

You can cancel the measurement process prematurely with any button; the pressure in the cuff is then quickly released.

Tip



Informing the patient

Go through the patient information sheet together with the patient before a 24-hour measurement. You can find the patient information sheet in A4 format at www.iem.de/mobil-o-graph/.

3.5 Positioning and behavior of the patient

Show the patient how the cuff must be positioned so that it is at the level of the right atrium during the measurement.

Inform the patient that they should adopt one of the positions shown in the table if possible when starting a blood pressure measurement.

1. Position	2. Position	3. Position
		

It is important that the patient:

- sits/stands/lies comfortably
- does **not** cross their legs
- places their feet flat on the floor (when sitting and standing)
- has their back and arms supported (when sitting and lying down)
- behaves calmly and does **not** speak

3.6 Technical data and ambient conditions

Measuring method	Oscillometric measuring method
Measuring pressure range:	Systolic 60 - 290 mmHg Diastolic 30 - 195 mmHg
Accuracy:	+/- 3 mmHg in the display range
Static pressure range:	0 - 300 mmHg
Pulse range:	30 - 240 beats per minute
Measuring intervals:	0, 1, 2, 3, 4, 5, 6, 10, 12, 15, 30 measurements per hour
Measurement protocols:	4 changeable interval groups
Storage capacity:	300 measurements (260 PWA)
Battery capacity:	> 300 measurements
Operating temperatures:	+10 °C - +40 °C

Operating air humidity:	15 % - 90 % relative humidity, non-condensing
Operating air pressure:	700 hPa - 1060 hPa
Transport/storage environment:	-25 °C - +70 °C, up to 90 % relative humidity, non-condensing
Dimensions:	128 x 75 x 30 mm
Weight:	approx. 180 g (without batteries)
Power supply:	2 Ni-MH rechargeable batteries with 1.2 V each and min. 1500mAh (AA, Mignon) 2 alkaline 1.5 V batteries (AA, Mignon)
Interfaces:	USB-C
Expected service life of the device	5 years
Expected service life of the cuff	6 months
Expected service life of the bag	5 years

Note

Ambient conditions:

- Extreme temperatures, humidity or air pressure can influence the measuring accuracy. Please observe the operating conditions.
- Extreme temperatures, humidity or altitude can affect the performance of the blood pressure monitor. Do not store the device near a fireplace or a radiant heater and do not expose it to extreme sunlight. Do not place the device next to a nebulizer or a kettle, as the condensation water can damage the device.
- The blood pressure monitor requires approx. 25 minutes to get from the minimum storage temperature of -25 °C to the operating temperature of +10 °C at an ambient temperature of +20 °C.
- The blood pressure monitor requires approx. 60 minutes to get from the maximum storage temperature of +70 °C to the operating temperature of +40 °C at an ambient temperature of +20 °C.

3.7 Symbols

Note

The symbols on the buttons are described in chapter 3.3.1 The buttons, p. 17.

Explanation of the nameplate symbols:

Symbol	Explanation
	Follow the instructions for use!
	The battery symbol indicates the type of recommended power supply.
	Manufacturer & Date of manufacture YYYY-MM-DD
	Type BF defibrillation-protected device
	The product should not be treated as normal household waste, but should be taken to a collection point for the recycling of electrical and electronic equipment. Further information can be obtained from your local authority, municipal waste disposal companies or specialist retailers.
	The device fulfills the basic requirements for labeling a medical device in accordance with Regulation (EU) 2017/745, certified by the Notified Body with the identification number XXXX (the digits after the CE mark depend on the Notified Body).
	Indicates that the product is a medical device
	MR-unsafe: The product presents hazards in MRI environments
	Serial number
	Protection class The blood pressure monitor is: **4**: Protected against solid foreign objects with a diameter ≥ 1.0 mm./Protected against access with a wire **2**: Protected against dripping water when the housing is tilted up to 15°.
	Metrological inspection label with next due date

4 Care and maintenance

Regular maintenance and care of the device is necessary to ensure that your Mobil-O-Graph® PWA continues to function perfectly.

The user (doctor) decides whether and when cleaning and/or disinfection of the blood pressure monitor, the cuff sleeve and the bag is necessary for hygienic reasons (e.g. after each use).

When the Mobil-O-Graph® PWA is not in use, remove the rechargeable batteries/batteries and observe the ambient conditions when storing the device, see 3.6 Technical data and ambient conditions, p. 23.

The device must be cleaned and disinfected before being handed over to the next patient.

NOTICE

Loss of warranty!

Opening the housing voids any warranty.

Do not open the device.

4.1 Cleaning

NOTICE

Damage to the device and accessories due to incorrect cleaning!

Incorrect, improper cleaning will render the device and/or the accessories unusable.

Clean the device and/or accessories as recommended by IEM GmbH.

NOTICE

Device damage due to penetrating liquid!

The device is not waterproof and will be rendered unusable due to penetrating liquid.

Do not immerse the cuff with bladder or the blood pressure monitor in disinfectant, water, or other liquids for cleaning.

However, if liquid has penetrated into the device or the device has been exposed to moisture, switch it off immediately. Remove the batteries and inform your healthcare provider without fail so that they can send it to the manufacturer for inspection.

Cleaning the blood pressure monitor

Only use a damp cotton cloth with lukewarm water to clean the blood pressure monitor.

Cleaning the cuff sleeve

- **Manual cleaning:**
 - Only use lukewarm water to clean the cuff sleeve and bag.
 - Do not use fabric softeners or other aids (e.g. hygiene rinses, textile deodorants). These agents can leave residues and damage the material!
- **Machine cleaning:**


Fig. 2: Removing the balloon from the cuff sleeve

a) Fully expand the cuff. If necessary, pull the cuff out of the hanger.

b) Remove the balloon with hose by pulling it out of the balloon compartment via the removal slot on the inside of the cuff.

- The cuff sleeve can be washed with a mild detergent in the washing machine at up to 30°C on a gentle spin cycle.
- The cuff sleeve is not suitable for tumble drying!

Tip

We also recommend the use of a laundry net.

The Velcro fastener should be closed before washing.

Do not wash the cuff sleeve with other products to avoid lint formation.

Cleaning balloon and hose manually

To clean the balloon and the hose, you must remove them from the cuff sleeve, Fig. 2: Removing the balloon from the cuff sleeve, p. 27.

Note

Clean the hose and balloon according to the cleaning instructions below:

- Use lukewarm water up to max. 30°C for cleaning.
 Never use harsh cleaning agents or cleaning agents containing solvents, as these can damage the surface of the blood pressure monitor and the cuff.
- Make sure that no water gets into the balloon or the hose!
- Never clean the balloon and hose mechanically!

1. Once the cuff, balloon and hose are completely dry, lay the cuff sleeve down and spread it out with the Velcro side facing down.
2. Slide the balloon over the removal slot into the balloon compartment and guide the hose through the small opening (from the inside of the cuff). Make sure that there are no creases when inserting!


Fig. 3: Inserting the balloon into the cuff sleeve

Note

- Make sure that the balloon size matches the cuff sleeve.
- The size information for the cuff sleeve can be found on the outside.

Cleaning the bag

- **Manual cleaning:**
 - Only use lukewarm water to clean the bag.
 - Do not use fabric softeners or other aids (e.g. hygiene rinses, textile deodorants). These agents can leave residues and damage the material!
- **Machine cleaning:**
 - The bag can be washed with a mild detergent in the washing machine at up to 30°C on a gentle spin cycle.
 - The bag is not suitable for tumble drying!

Tip

We also recommend the use of a laundry net.

The Velcro fastener should be closed before washing.

Do not wash the bag with other products to avoid lint formation.

4.2 Disinfection

IEM GmbH has tested the following surface disinfectants for the disinfection of the blood pressure monitor and the approved accessories:

- Meliseptol ® New Formular (B. Braun)
- Mikrozid ® universal liquid (Schülke)
- Mikrozid ® AF liquid (Schülke)

When using other disinfectants that have not been tested by IEM, the user is responsible for providing proof that they are harmless to use.

Never use disinfectants that leave residues on the product or that are not suitable for contact with the skin.

Tip

For further information on the surface disinfectants, their properties, compositions and efficacy spectra, please refer to the respective manufacturer's information!

Disinfecting blood pressure monitor and approved accessories

1. Spray the blood pressure monitor and accessories with one of the disinfectants listed above. Make sure that no liquid penetrates the blood pressure monitor!
2. Always allow the product to dry completely.

NOTICE

Device damage due to penetrating liquid!

The device is not waterproof and will be rendered unusable due to penetrating liquid.

Do not immerse the cuff with bladder or the blood pressure monitor in disinfectant, water, or other liquids for cleaning.

NOTICE

However, if liquid has penetrated into the device or the device has been exposed to moisture, switch it off immediately. Remove the batteries and inform your healthcare provider without fail so that they can send it to the manufacturer for inspection.

4.3 Maintenance schedule

Every 2 years:

As proof of continuous compliance to “Basic Requirements” pursuant to the EU regulation 2017/745, the Mobil-O-Graph® PWA blood pressure monitor must undergo metrological checks every two years. In certain countries, this requirement may be regulated by national laws or regulations.

We offer to take over the metrological inspection and the complete service.

This service includes:

- Metrological inspection
- Software update (if available)
- Function test: Electronics, pump and pneumatic circuit

Apart from the metrological inspection, no further maintenance measures are required for electromagnetic compatibility. The due date of the next metrological inspection can be found on the metrological inspection label.

5 Troubleshooting

Note

In the event of an error measurement, the device starts a new measurement after 3 minutes, except for the initial measurement.

5.1 Main sources of error

The following can cause error measurements or unwanted results:

- The patient’s arm moving during measurement
- Device being switched off (e.g. at night)
- Incorrect cuff size put on
- The cuff slipping while being worn
- Lack of successful manual measurement in the doctor's practice
- Check that the connection cable is connected correctly and that **co** appears on the display
- Not fully charged, incorrectly charged or outdated rechargeable batteries used
- Kinking or knotting of the cuff tube **(during the sleep phase)**
- External interference factors such as, for example, physical activity or driving or using public transport during measurement may result in so-called motion artefacts or incorrect measurements.

5.2 Error messages

Display	Troubleshooting
bAtt E005	Replace the batteries.
-SToP- E008	*Measurement interrupted by pressing a button*
rEboot	If rEboot appears on the display but the device then behaves normally, you do not need to do anything. If this message precedes a device error: Send the device immediately to your specialist dealer or directly to IEM GmbH for checking.

Display	Troubleshooting
E001	Internal error: Send the device to your specialist dealer or directly to IEM GmbH for checking.
E002	The arm was moved during measurement. The arm must be kept as still as possible during the measurement. The patient has arrhythmia. Blood pressure outside measurement range. If this message continues, the device is not suitable for the patient.
E003	Blood pressure outside measurement range. If this message continues, the device is not suitable for the patient. The arm was moved during measurement. The arm must be kept as still as possible during the measurement. Pneumatic problems. If the problem occurs frequently, send the device to your specialist dealer or directly to IEM GmbH for checking.
E004	The data transfer cable is not plugged in correctly. Check the USB-C connector for damage.
E006	The cuff is not properly connected or is affected by a kink or blockage. Pneumatic problems. If the problem occurs frequently, send the device to your specialist dealer or directly to IEM GmbH for checking.
E007	The device's memory is full (maximum 300 entries).Delete the data from the device. Ensure in advance that the existing data has been saved and then reinitialize the device.
E009	Pneumatic problems. If the problem occurs frequently, send the device to your specialist dealer or directly to IEM GmbH for checking.
E010	These error messages all indicate a serious error.Send the device to your specialist dealer or directly to IEM GmbH for checking.

6 Rechargeable batteries/batteries

Important

To ensure safe operation of the device, observe and follow the relevant safety and hazard information: Batteries/rechargeable batteries, p. 14

Note

Inserting/removing the batteries/rechargeable batteries
Switch the device off before inserting or removing the batteries/rechargeable batteries.

6.1 Operation with rechargeable batteries

IEM GmbH recommends high-quality "ready to use" NiMH rechargeable batteries with at least 1500 mAh, which have a very low self-discharge rate. We recommend using only such batteries!

Note

Batteries are subject to ageing.
Batteries that are damaged or with which a 24-hour measurement is no longer possible must be replaced immediately.

Note

Do not use rechargeable batteries that have been stored above 45 °C or below 0 °C.
If the operating time of the rechargeable batteries drops significantly, replace all rechargeable batteries at the same time.
Never use old used rechargeable batteries together with new unused rechargeable batteries!

Note

Observe minimum voltage
Always ensure that the batteries/rechargeable batteries have a sufficient minimum voltage of **at least 2.4 V** to guarantee operation of the blood pressure monitor.
Always use batteries/rechargeable batteries with full capacity for a 24-hour blood pressure measurement (minimum voltage required for a 24-hour measurement series: 2.6 V).

Tips:

- Fully charge the rechargeable batteries before using them for the first time.
- Please note that NiMH rechargeable batteries only reach their full charging capacity after the 4th charging cycle.
- Recharge the batteries if they have not been used for a long time.
- To protect your rechargeable batteries, avoid deep discharging.

Tip

When charging rechargeable batteries, follow the instructions for use provided by the manufacturer of the used charger.

6.2 Operation with ALKALINE batteries

The Mobil-O-Graph® PWA can also be operated with alkaline manganese batteries. These batteries are not rechargeable. When fully charged, two alkaline manganese batteries have a higher charging voltage (at least 3.1 V) compared to two NiMH batteries (at least 2.6 V). The battery indicator in the blood pressure monitor will therefore also show a higher voltage.

7 Disposal

Blood Pressure Monitor

The symbol on the product or the packaging indicates that this product must not be treated like normal household waste.

In accordance with the Elektro- und Elektronikgerätegesetz (ElektroG, [German Electronic and Electrical Equipment Act]), IEM GmbH will take back old electrical devices at any time. They will then either be reused or disposed of using a recycling company in accordance with statutory regulations.

The IEM WEEE reg. no. is: **DE 85825018**

The customer themselves must bear costs incurred that are directly related to the return of old electrical devices, such as transport costs, import duties, customs clearance costs and any other additional charges related to logistics. The customer must organize the return of old electrical devices to IEM themselves.

Please use the following address for returning old electrical devices:

IEM GmbH
Gewerbepark Brand 42
52078 Aachen
Germany

Approved accessories

Approved, non-electronic accessories (e.g. IEM ABPM cuff and device bag) can be disposed of in normal household waste.

Data Protection during Disposal

IEM informs all end users of old electronic and electrical devices that they are personally responsible for deleting personal data on the old devices to be disposed of.

Batteries/Rechargeable Batteries

Batteries do not belong in household waste. As a consumer, you are legally obliged to return used batteries. You can drop off old batteries at the public collection points in your community or anywhere where batteries of the relevant type are sold.

	Li	Battery contains lithium
	Al	Battery contains alkali
	Mn	Battery contains manganese

8 Warranty and repair conditions

Warranty:

IEM GmbH provides a two-year warranty on the hardware, i.e. the blood pressure monitor itself. The warranty is no longer valid if the device has been opened, misused or negligently or intentionally damaged. Likewise, the warranty is no longer valid if the device has been repaired by an unauthorised body. Only the NiMH rechargeable batteries can be replaced by the operator.

The normal wear and tear of rechargeable batteries, transmission cables or cuffs including tube cannot be claimed under the warranty. In the scope of the warranty, IEM is only responsible for defects that already existed when the product was handed over to the customer (and possibly only showed up later).

Furthermore, the warranty does not apply if the defect is due to violation of the instructions for use, improper handling, moisture, humidity or extreme heat or climate conditions or short-term fluctuations in corresponding influences or is due to corrosion, oxidation, unauthorised intervention or connection attempts, unauthorised opening or repair, attempts at repair with non- approved replacement parts, incorrect operation, improper installation, accidents, forces of nature, spillage of food or drink, chemical effects or other external influences which IEM cannot influence (among other things, defects in consumables such as batteries and rechargeable batteries, for example, which inevitably only have a limited operating life, are excluded), unless the defect is directly due to a material defect, design fault or manufacturing error.

Repairs:

If the device malfunctions or if you experience continuous measurement errors, contact your Mobil-O-Graph® PWA specialist supplier or IEM GmbH directly for repairs and shipping instructions.

Note

Do not open the housing.

Once the device has been opened, any warranty will no longer be valid.

Metrological control:

As proof of continuous compliance to “Basic Requirements” pursuant to the EU regulation 2017/745, the Mobil-O-Graph® PWA blood pressure monitor must undergo metrological checks every two years. In certain countries, this requirement may be regulated by national laws or regulations.

Liability clause:

In all cases where IEM is obliged to pay damages or reimburse expenses on the basis of contractual or statutory claims, IEM is liable only to the extent that its executives and vicarious agents are guilty of intent or gross negligence. Strict liability according to the Product Liability Act remains unaffected. Liability for the culpable breach of essential contractual obligations also remains unaffected, however, except in the cases of clauses 1 and 2, liability is limited to foreseeable damage which is typical for this type of contract.

A change in the burden of proof to the disadvantage of the customer is not connected with the above regulations. Liability for quality or durability guarantees or fraudulent concealment of defects is not affected by the provisions outlined above.

Appendix

A Appendix 1: EMC Guidelines and Manufacturer's Declaration

Guidelines and Manufacturer's Declaration – electromagnetic emissions		
The blood pressure measuring device Mobil-O-Graph® PWA is intended for use in the environment specified below. The customer or user of the blood pressure measuring device Mobil-O-Graph® PWA should ensure that it is used in such an environment.		
Emitted interference measurement	**Compliance**	**Electromagnetic environment – guidelines**
RF interference emissions according to CISPR 11	Group 1	The blood pressure measuring device Mobil-O-Graph® PWA only uses RF power for its internal functioning. Its RF emission is therefore very low and it is unlikely that any neighboring electronic device will experience any interference.
RF interference emissions according to CISPR 11	Class B	The blood pressure measuring device Mobil-O-Graph® PWA is suitable for use in residential areas as well as in other facilities that are directly connected to a public supply network that also supplies buildings used for residential purposes.
RF interference emissions according to CISPR 25	Not applicable	
Harmonic oscillations according to IEC 61000-3-2	Not applicable	
Voltage fluctuations/flickers according to IEC 61000-3-3	Not applicable	

Guidelines and manufacturer's declaration – electromagnetic immunity			
The blood pressure measuring device Mobil-O-Graph® PWA is intended for use in the environment specified below. The customer or user of the blood pressure measuring device Mobil-O-Graph® PWA should ensure that it is used in such an environment.			
Immunity test	**Test level**	**Compliance level**	**Electromagnetic environment – guidelines**
Electrostatic discharge (ESD) according to IEC 61000-4-2	+- 8 kV contact discharge +- 15 kV air discharge	+- 8 kV contact discharge +- 15 kV air discharge	Floors should consist of wood or cement or be covered with ceramic tiles. If the floor consists of synthetic materials, relative humidity must be at least 30 %.
Fast transient electrical disturbance/bursts according to IEC 61000-4-4	± 1 kV 100 kHz repetition rate	± 1 kV 100 kHz repetition rate	
Surge voltages (surges) according to IEC 61000-4-5		Not applicable	The Mobil-O-Graph® PWA does not have an AC power supply

Guidelines and manufacturer's declaration – electromagnetic immunity			
Magnetic field at supply frequency (50/60 Hz) according to IEC 61000-4-8	30 A/m	30 A/m	Magnetic fields at mains frequency should match the typical values found in business or hospital environments.
Voltage dips, short interruptions and fluctuations in supply voltage according to IEC 61000-4-11		Not applicable	The Mobil-O-Graph® PWA does not have an AC power supply
Radiated RF disturbance according to IEC 61000-4-3	10 V/m 80 MHz to 2,7 GHz	10 V/m	
Conducted disturbances according to IEC 61000-4-6		Not applicable	
Radiated RF disturbance according to IEC 61000-4-3	380 - 390 MHz 27 V/m; PM 50%; 18 Hz	380 - 390 MHz 27 V/m; PM 50%; 18 Hz	
	430 - 470 MHz 28 V/m; (FM ±5 kHz, 1 kHz sine) PM; 18 Hz	430 - 470 MHz 28 V/m; (FM ±5 kHz, 1 kHz sine) PM; 18 Hz	
	704 - 787 MHz 9 V/m; PM 50%; 217 Hz	704 - 787 MHz 9 V/m; PM 50%; 217 Hz	
	800 - 960 MHz 28 V/m; PM 50%; 18 Hz	800 - 960 MHz 28 V/m; PM 50%; 18 Hz	
	1700 - 1990 MHz 28 V/m; PM 50%; 217 Hz	1700 - 1990 MHz 28 V/m; PM 50%; 217 Hz	
	2400 - 2570 MHz 28 V/m; PM 50%; 217 Hz	2400 - 2570 MHz 28 V/m; PM 50%; 217 Hz	
	5100 - 5800 MHz 9 V/m; PM 50%; 217 Hz	5100 - 5800 MHz 9 V/m; PM 50%; 217 Hz	
Magnetic near fields according to IEC 61000-4-39	30 kHz, modulation, CW 8 A/m 134.2 kHz, 2.1 kHz, PM13 65 A/m 13.56 MHz, 50 kHz, PM 7.5 A/m	30 kHz, modulation, CW 8 A/m 134.2 kHz, 2.1 kHz, PM13 65 A/m 13.56 MHz, 50 kHz, PM 7.5 A/m	

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